

EF-17
Ortho-SUV MINI FRAME
(Computer-Assisted)
EXTERNAL FIXATOR
SET

Ortho-SUV MINI FRAME FIXATOR

Ortho SUV Mini is a new generation computer-assisted external fixation device, a hexapod system for treatment of pediatric deformities. Ortho-SUV Mini Frame is an advanced version of the traditional Ring Fixator and developed as a single-step solution for treating complex deformities and fractures.

All the unique advantages of "elder brother", OSF-7, are kept:


- *Adjustable and useable with all kind of rings
- *Offline image-based software utilization with X-ray loaded calculations and error checking
- *12 parameters are also manually measured
- *Single-step correction for Multi Planar Deformity & Fracture Reduction
- *Non-predetermined places for Struts Fixation and non-predetermined directions for the rings
- *Dual Mode Struts – Fast & Gradual Correction Mode
- *Possibility to use two structures at risk (SAR): bone- and soft tissue depended.

Indications: Diaphyseal and metadiaphyseal long bone deformities of any complicity. Complex deformities of mid-foot and hindfoot. Rigid contractures and old dislocations and subluxations of the knee, elbow, wrist, and ankle joints. Diaphyseal and metadiaphyseal long bone fractures.



Implants & Instruments

TST Code	Barcode	Description	Qty
00000506100	8680858451906	ORTHO-SUV 1. TRAY	1
52450010075	8680858451487	MINI STRUT LARGE 75 MM	6
52450010100	8680858451494	MINI STRUT LARGE 100 MM	6
52450010125	8680858451500	MINI STRUT LARGE 125 MM	6
52450010150	8680858451517	MINI STRUT LARGE 150 MM	6
52450010200	8680858451524	MINI STRUT LARGE 200 MM	6
52450010250	8680858451531	MINI STRUT LARGE 250 MM	6
52450020075	8680858451548	MINI STRUT MEDIUM 75 MM	6
52450020100	8680858451555	MINI STRUT MEDIUM 100 MM	6
52450020125	8680858451562	MINI STRUT MEDIUM 125 MM	6
52450020150	8680858451579	MINI STRUT MEDIUM 150 MM	6
52450020200	8680858451586	MINI STRUT MEDIUM 200 MM	6
52450020250	8680858451593	MINI STRUT MEDIUM 250 MM	6
52450030075	8680858451609	MINI STRUT SMALL 75 MM	6
52450030100	8680858451616	MINI STRUT SMALL 100 MM	6
52450030125	8680858451623	MINI STRUT SMALL 125 MM	6
52450030150	8680858451630	MINI STRUT SMALL 150 MM	6
52450030200	8680858451647	MINI STRUT SMALL 200 MM	6
52450030250	8680858451654	MINI STRUT SMALL 250 MM	6



Implants & Instruments

No	TST Code	Barcode	Description	Qty
1	00000506200	8680858451913	ORTHO-SUV 2. TRAY	1
2	52870000001	8680858444526	X-RAY MARKER OF STRUT NO.1	1
	52870000002	8680858444533	X-RAY MARKER OF STRUT NO.2	1
	52870000003	8680858444540	X-RAY MARKER OF STRUT NO.3	1
	52870000004	8680858444557	X-RAY MARKER OF STRUT NO.4	1
	52870000005	8680858444564	X-RAY MARKER OF STRUT NO.5	1
	52870000006	8680858444571	X-RAY MARKER OF STRUT NO.6	1
3	52920000001	8680858444588	LABEL OF STRUT NO.1	1
	52920000002	8680858444595	LABEL OF STRUT NO.2	1
	52920000003	8680858444601	LABEL OF STRUT NO.3	1
	52920000004	8680858444618	LABEL OF STRUT NO.4	1
	52920000005	8680858444625	LABEL OF STRUT NO.5	1
	52920000006	8680858444632	LABEL OF STRUT NO.6	1
4	53010000800	8680858444649	SRC STRAIGHT PLATE	6
5	53010000900	8680858444656	SRC Z PLATE	3
6	04050000008	8699931028379	WRENCH Ø 8 MM	2
7	08044000012	8699931015744	WRENCH Ø 12 MM	2
8	04015000040	8698673490055	ALLEN WRENCH Ø 4 MM	1
9	05241000001	8680858444663	TRIANGULAR MEASUREMENT DEVICE	1
10	05240000300	8680858444670	RULER 300 MM	1
	00554273135	8680858440986	CONTAINER 554X273X135 MM AL.	1

