

EF-08
ULTRA MODULAR EXTERNAL FIXATOR
(UMEF) SET



UMEF FIXATOR


UMEF “Ultra Modular External Fixator” is a unique system that allows for a rigid osteosynthesis application, developed aiming at providing an effective continuation of the stabilization in the period of time to fracture healing after performing the fracture repositioning. The system consists of rods, clamps, and pins. Rod-rod clamps serve to connect two rods to each other and rod – pin clamps that provide to connect between rod and pin. UMEF is lightweight owing to its carbon fibre structure of rods and aluminium made clamps.

Indications: For fracture fixation of lower and upper limbs.

rod-pin clamp double

rod-rod clamp


rod-pin clamp single



Implants & Instruments

No	TST Code	Barcode	Description	Qty
1	00000701100	8699931017151	UMEF HUMERUS/ULNA RADIUS 1. TRAY	1
2	00000601300	8699931017182	3-4 MM SCHANZ BOX	1
	20910602003	8698673455887	DAF SCHANZ PIN CORTICAL THREADED 3X20X60	8
	20910802003	8698673455894	DAF SCHANZ PIN CORTICAL THREADED 3X20X80	8
	20910002503	8698673455917	DAF SCHANZ PIN CORTICAL THREADED 3X25X100	12
	20910202003	8698673455900	DAF SCHANZ PIN CORTICAL THREADED 3X20X120	8
	20910702004	8698673458604	DAF SCHANZ PIN CORTICAL THREADED 4X20X70	4
	20910802004	8698673455962	DAF SCHANZ PIN CORTICAL THREADED 4X20X80	8
	20910002004	8698673473317	DAF SCHANZ PIN CORTICAL THREADED 4X20X100	8
	20910202504	8698673456020	DAF SCHANZ PIN CORTICAL THREADED 4X25X120	8
	20910142504	8698673456037	DAF SCHANZ PIN CORTICAL THREADED 4X25X140	8
3	04051300004	8699931010107	UMEF SCHANZ GUIDE 4 MM	1
4	04051302004	8699931010114	UMEF SCHANZ GUIDE DOUBLE 4 MM	1
5	04015000300	8698673492493	3 MM SCHANZ PIN T DRIVER	1
6	04015000400	8698673492509	4 MM SCHANZ PIN T DRIVER	1

* TST Code	Barcode	Description
20920062003	8698673455306	SCHANZ PIN CORTICAL THREADED 3X20X60 MM
20920082003	8698673455313	SCHANZ PIN CORTICAL THREADED 3X20X80 MM
20920102503	8698673455320	SCHANZ PIN CORTICAL THREADED 3X25X100 MM
20921002003	8699931022018	SCHANZ PIN CORTICAL THREADED 3X20X120 MM
20920072004	8698673458741	SCHANZ PIN CORTICAL THREADED 4X20X70 MM
20920082004	8698673455337	SCHANZ PIN CORTICAL THREADED 4X20X80 MM
20920102004	8698673455344	SCHANZ PIN CORTICAL THREADED 4X20X100 MM
20920122504	8698673455368	SCHANZ PIN CORTICAL THREADED 4X25X120 MM
20921402504	8699931028607	SCHANZ PIN CORTICAL THREADED 4X25X140 MM



Implants & Instruments

No	TST Code	Barcode	Description	Qty
7	00000701200	8699931017168	UMEF HUMERUS/ULNA RADIUS 2.TRAY	1
8	40523300001	8698673468382	UMEF HUMERUS ROD-PIN CLAMP SINGLE	5
9	40523300002	8698673468375	UMEF HUMERUS ROD-PIN CLAMP DOUBLE	5
10	40523200001	8698673468399	UMEF HUMERUS ROD-ROD CLAMP	5
11	04050000009	8699931018028	WRENCH 9 MM	2
12	00000701300	8699931017175	UMEF HUMERUS 3.TRAY	1
13	40563108060	8698673468313	UMEF HUMERUS ROD FIBER CARBON 8X60 MM	5
14	40563108120	8698673468276	UMEF HUMERUS ROD FIBER CARBON 8X120 MM	5
15	40563108180	8699931004557	UMEF HUMERUS ROD FIBER CARBON 8X180 MM	5
16	40563108240	8698673468290	UMEF HUMERUS ROD FIBER CARBON 8X240 MM	5
17	40524000001	8698673468412	UMEF ULNA RADIUS COMPR. DISTR.UNIT	1
18	40520300001	8699931018974	UMEF ULNA RADIUS COMPR. DISTR.UNIT(SHORT)	1
	00500160120	8699931017298	500X160X120 MM CONTAINER	1

