

TRAUMA TIBIA & FIBULA



Distal Posterior Lateral Tibial locking plate

Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
11.707.25.405	L	5	62.2
11.707.25.505	R	5	62.2
11.707.25.407	L	7	90.2
11.707.25.507	R	7	90.2



TS Distal Tibial Locking Plate (Anterior Column)

Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
11.707.25.605	L	5	78
11.707.25.705	R	5	78
11.707.25.607	L	7	106
11.707.25.707	R	7	106



Locking Compression Plate III (T-shaped)

Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
15.307.26.204	L	4	97.5
15.307.26.205	L	5	110.5
15.307.26.206	L	6	123.5
15.307.26.208	L	8	149.5
15.307.26.210	L	10	175.5
15.307.26.212	L	12	201.5
15.307.26.304	R	4	97.5
15.307.26.305	R	5	110.5
15.307.26.306	R	6	123.5
15.307.26.308	R	8	149.5
15.307.26.310	R	10	175.5
15.307.26.312	R	12	201.5



Locking Compression Plate III (L-shaped)

Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
15.307.27.404	L	4	81.5
15.307.27.405	L	5	94.5
15.307.27.406	L	6	107.5
15.307.27.408	L	8	133.5
15.307.27.410	L	10	159.5
15.307.27.412	L	12	185.5
15.307.27.504	R	4	81.5
15.307.27.505	R	5	94.5
15.307.27.506	R	6	107.5
15.307.27.508	R	8	133.5
15.307.27.510	R	10	159.5
15.307.27.512	R	12	185.5



Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
15.307.24.004	L	4	116
15.307.24.006	L	6	142
15.307.24.008	L	8	168
15.307.24.010	L	10	194
15.307.24.012	L	12	220
15.307.24.014	L	14	246
15.307.24.104	R	4	116
15.307.24.106	R	6	142
15.307.24.108	R	8	168
15.307.24.110	R	10	194
15.307.24.112	R	12	220
15.307.24.114	R	14	246



Distal Lateral Tibial Locking Compression Plate

Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
11.107.37.005	L	5	81
11.107.37.007	L	7	107
11.107.37.009	L	9	133
11.107.37.011	L	11	159
11.107.37.013	L	13	185
11.107.37.105	R	5	81
11.107.37.107	R	7	107
11.107.37.109	R	9	133
11.107.37.111	R	11	159
11.107.37.113	R	13	185



Anatomical Dital Lateral Fibular Locking Compression Plate

Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
15.307.14.003	L	3	73
15.307.14.004	L	4	86
15.307.14.005	L	5	99
15.307.14.006	L	6	112
15.307.14.007	L	7	125
15.307.14.009	L	9	151
15.307.14.011	L	11	177
15.307.14.013	L	13	203
15.307.14.103	R	3	73
15.307.14.104	R	4	86
15.307.14.105	R	5	99
15.307.14.106	R	6	112
15.307.14.107	R	7	125
15.307.14.109	R	9	151
15.307.14.111	R	11	177
15.307.14.113	R	13	203



Distal Lateral Fibular Locking Compression Plate

Ti:

PRODUCT CODE	L/R	HOLES	LENGTH
11.107.15.004	L	4	92
11.107.15.006	L	6	116
11.107.15.008	L	8	140
11.107.15.010	L	10	164
11.107.15.012	L	12	188
11.107.15.014	L	14	212
11.107.15.104	R	4	92
11.107.15.106	R	6	116
11.107.15.107	R	7	128
11.107.15.108	R	8	140
11.107.15.110	R	10	164
11.107.15.112	R	12	188
11.107.15.114	R	14	212



Distal Posterior Lateral Fibular Locking Compression Plate

Ti:

PRODUCT CODE	DIAMETER	LENGTH
15.307.40.003	3	84
15.307.40.005	5	110
15.307.40.007	7	136
15.307.40.009	9	162
15.307.40.011	11	188
15.307.40.013	13	214



Proximal Lateral Tibial Locking Compression Plate III

Ti:

PRODUCT CODE	L/R	DIAMETER	LENGTH
11.107.39.403	L	3	76
11.107.39.404	L	4	108
11.107.39.406	L	6	140
11.107.39.408	L	8	172
11.107.39.410	L	10	204
11.107.39.412	L	12	236
11.107.39.503	R	3	76
11.107.39.504	R	4	108
11.107.39.506	R	6	140
11.107.39.508	R	8	172
11.107.39.510	R	10	204
11.107.39.512	R	12	236

(NOT IN SA)



Tibia Locking Compression Plate

Ti:

PRODUCT CODE	DIAMETER	LENGTH
15.307.02.006	6	109
15.307.02.007	7	127
15.307.02.008	8	145
15.307.02.009	9	163
15.307.02.010	10	181
15.307.02.012	12	217
15.307.02.014	14	253
15.307.02.016	16	289



Locking Compression Plate IV (T-shaped)

Ti:

PRODUCT CODE	L/R	DIAMETER	LENGTH
15.307.26.404	L	4	105
15.307.26.405	L	5	123
15.307.26.406	L	6	141
15.307.26.408	L	8	177
15.307.26.410	L	10	213
15.307.26.412	L	12	249
15.307.26.504	R	4	105
15.307.26.505	R	5	123
15.307.26.506	R	6	141
15.307.26.508	R	8	177
15.307.26.510	R	10	213
15.307.26.512	R	12	249



Proximal Medial Tibial Osteotomy Locking Compression Plate(T-Shape)

(NOT IN SA)

Ti:

PRODUCT CODE	L/R	DIAMETER	LENGTH
11.109.33.004	L	4	115
11.109.33.104	R	4	115



HA3.5mm Locking Screw

- Full Thread
- Self Tapping
-  T15



Ti:

PRODUCT CODE	LENGTH
11.107.51.110	10
11.107.51.112	12
11.107.51.114	14
11.107.51.116	16
11.107.51.118	18
11.107.51.120	20
11.107.51.122	22
11.107.51.124	24
11.107.51.126	26
11.107.51.128	28

PRODUCT CODE	LENGTH
11.107.51.130	30
11.107.51.132	32
11.107.51.134	34
11.107.51.135	35
11.107.51.136	36
11.107.51.138	38
11.107.51.140	40
11.107.51.142	42
11.107.51.144	44
11.107.51.145	45

PRODUCT CODE	LENGTH
11.107.51.146	46
11.107.51.148	48
11.107.51.150	50
11.107.51.155	55
11.107.51.160	60
11.107.51.165	65
11.107.51.170	70
11.107.51.175	75
11.107.51.180	80
11.107.51.185	85

HA3.5mm Cortex Screw

- Full Thread
- Self Tapping
-  SW2.5



Ti:

PRODUCT CODE	LENGTH
11.105.03.010	10
11.105.03.012	12
11.105.03.014	14
11.105.03.016	16
11.105.03.018	18
11.105.03.020	20
11.105.03.022	22
11.105.03.024	24

PRODUCT CODE	LENGTH
11.105.03.026	26
11.105.03.028	28
11.105.03.030	30
11.105.03.032	32
11.105.03.034	34
11.105.03.036	36
11.105.03.038	38
11.105.03.040	40

PRODUCT CODE	LENGTH
11.105.03.042	42
11.105.03.044	44
11.105.03.046	46
11.105.03.048	48
11.105.03.050	50
11.105.03.055	55
11.105.03.060	60

HA4.5mm Cortex Screw

- Full Thread
- Self Tapping
-  SW3.5



Ti:

PRODUCT CODE	LENGTH
11.105.08.016	16
11.105.08.018	18
11.105.08.020	20
11.105.08.022	22
11.105.08.024	24
11.105.08.026	26
11.105.08.028	28
11.105.08.030	30

PRODUCT CODE	LENGTH
11.105.08.032	32
11.105.08.034	34
11.105.08.036	36
11.105.08.038	38
11.105.08.040	40
11.105.08.042	42
11.105.08.044	44
11.105.08.046	46

PRODUCT CODE	LENGTH
11.105.08.048	48
11.105.08.050	50
11.105.08.052	52
11.105.08.054	54
11.105.08.056	56
11.105.08.058	58
11.105.08.060	60

HA5.0mm Locking Screw

- Full Thread
- Self Tapping
-  T25



Ti:

PRODUCT CODE	LENGTH
11.107.55.114	14
11.107.55.116	16
11.107.55.118	18
11.107.55.120	20
11.107.55.122	22
11.107.55.124	24
11.107.55.126	26
11.107.55.128	28
11.107.55.130	30
11.107.55.132	32
11.107.55.134	34


PRODUCT CODE	LENGTH
11.107.55.136	36
11.107.55.138	38
11.107.55.140	40
11.107.55.142	42
11.107.55.144	44
11.107.55.146	46
11.107.55.148	48
11.107.55.150	50
11.107.55.155	55
11.107.55.160	60
11.107.55.165	65

PRODUCT CODE	LENGTH
11.107.55.170	70
11.107.57.075	75
11.107.57.080	80
11.107.57.085	85
11.107.57.090	90
11.107.57.095	95
11.107.57.100	100
11.107.57.105	105
11.107.57.110	110

HB4.0mm Cancellous Screw

- Full Thread

-  SW2.5



Ti:

PRODUCT CODE	LENGTH
11.105.20.012	12
11.105.20.014	14
11.105.20.016	16
11.105.20.018	18
11.105.20.020	20
11.105.20.022	22
11.105.20.024	24
11.105.20.026	26
11.105.20.028	28

PRODUCT CODE	LENGTH
11.105.20.030	30
11.105.20.032	32
11.105.20.034	34
11.105.20.036	36
11.105.20.038	38
11.105.20.040	40
11.105.20.042	42
11.105.20.044	44
11.105.20.046	46

PRODUCT CODE	LENGTH
11.105.20.048	48
11.105.20.050	50
11.105.20.052	52
11.105.20.054	54
11.105.20.056	56
11.105.20.058	58
11.105.20.060	60

HB4.0mm Cancellous Screw

- Half Thread

-  SW2.5



Ti:

PRODUCT CODE	LENGTH
11.105.21.014	14
11.105.21.016	16
11.105.21.018	18
11.105.21.020	20
11.105.21.022	22
11.105.21.024	24
11.105.21.026	26
11.105.21.028	28

PRODUCT CODE	LENGTH
11.105.21.030	30
11.105.21.032	32
11.105.21.034	34
11.105.21.036	36
11.105.21.038	38
11.105.21.040	40
11.105.21.042	42
11.105.21.044	44

PRODUCT CODE	LENGTH
11.105.21.046	46
11.105.21.048	48
11.105.21.050	50
11.105.21.052	52
11.105.21.054	54
11.105.21.056	56
11.105.21.058	58
11.105.21.060	60

HB6.5mm Cancellous Screw

- Full Thread

-  SW3.5



Ti:

PRODUCT CODE	LENGTH
11.105.24.025	25
11.105.24.030	30
11.105.24.035	35
11.105.24.040	40
11.105.24.045	45
11.105.24.050	50
11.105.24.055	55

PRODUCT CODE	LENGTH
11.105.24.060	60
11.105.24.065	65
11.105.24.070	70
11.105.24.075	75
11.105.24.080	80
11.105.24.085	85
11.105.24.090	90

PRODUCT CODE	LENGTH
11.105.24.095	95
11.105.24.100	100
11.105.24.105	105
11.105.24.110	110
11.105.24.115	115
11.105.24.120	120

HB6.5mm Cancellous Screw

- Half Thread

-  SW3.5



Ti:

PRODUCT CODE	THREAD	TOTAL
11.105.26.030	16	30
11.105.26.035	16	35
11.105.26.040	16	40
11.105.26.045	16	45
11.105.26.050	16	50
11.105.26.055	16	55
11.105.26.060	16	60
11.105.26.065	16	65
11.105.26.070	16	70

PRODUCT CODE	THREAD	TOTAL
11.105.26.075	16	75
11.105.26.080	16	80
11.105.26.085	16	85
11.105.26.090	16	90
11.105.26.095	16	95
11.105.26.100	16	100
11.105.26.105	16	105
11.105.26.110	16	110

ANKLE LOCKING PLATE

Hollow Tibial Osteotomy Plate

Ti:

PRODUCT CODE	L/R	HOLES
11.108.13.009	L	9
11.108.13.010	L	10
11.108.13.109	R	9
11.108.13.110	R	10



Anterior Ankle Fusion Plate

Ti:

PRODUCT CODE	L/R	HOLES
11.108.14.003	L	3
11.108.14.004	L	4
11.108.14.103	R	3
11.108.14.104	R	4



HA3.5mm Locking Screw

-Full Thread

-Self Tapping

-  T15



Ti:

PRODUCT CODE	LENGTH
11.108.54.020	20
11.108.54.022	22
11.108.54.024	24
11.108.54.026	26
11.108.54.028	28
11.108.54.030	30

PRODUCT CODE	LENGTH
11.108.54.032	32
11.108.54.034	34
11.108.54.036	36
11.108.54.038	38
11.108.54.040	40
11.108.54.042	42

PRODUCT CODE	LENGTH
11.108.54.044	44
11.108.54.046	46
11.108.54.048	48
11.108.54.050	50

HA3.5mm Cortex Screw

-Full Thread

-Self Tapping

-  T15



Ti:

PRODUCT CODE	LENGTH
11.108.55.020	20
11.108.55.022	22
11.108.55.024	24
11.108.55.026	26
11.108.55.028	28
11.108.55.030	30

PRODUCT CODE	LENGTH
11.108.55.032	32
11.108.55.034	34
11.108.55.036	36
11.108.55.038	38
11.108.55.040	40
11.108.55.042	42

PRODUCT CODE	LENGTH
11.108.55.044	44
11.108.55.046	46
11.108.55.048	48
11.108.55.050	50

X-perts Tibia Intramedullary Nail



Ti:

PRODUCT CODE	DIAMETER	LENGTH
11.106.17.255	8.5	255
11.106.17.270	8.5	270
11.106.17.285	8.5	285
11.106.17.300	8.5	300
11.106.17.315	8.5	315
11.106.17.330	8.5	330
11.106.17.345	8.5	345
11.106.17.360	8.5	360
11.106.17.375	8.5	375
11.106.18.255	9	255
11.106.18.270	9	270
11.106.18.285	9	285
11.106.18.300	9	300
11.106.18.315	9	315
11.106.18.330	9	330
11.106.18.345	9	345
11.106.18.360	9	360
11.106.18.375	9	375

11.106.19.255	10	255
11.106.19.270	10	270
11.106.19.285	10	285
11.106.19.300	10	300
11.106.19.315	10	315
11.106.19.330	10	330
11.106.19.345	10	345
11.106.19.360	10	360
11.106.19.375	10	375
11.106.20.255	11	255
11.106.20.270	11	270
11.106.20.285	11	285
11.106.20.300	11	300
11.106.20.315	11	315
11.106.20.330	11	330
11.106.20.345	11	345
11.106.20.360	11	360
11.106.20.375	11	375

Tibia Locking Screw

- Full Thread
- Self Tapping
-  SW4.5



Ti:

PRODUCT CODE	LENGTH
11.106.26.024	24
11.106.26.026	26
11.106.26.028	28
11.106.26.030	30
11.106.26.032	32
11.106.26.034	34
11.106.26.036	36
11.106.26.038	38
11.106.26.040	40
11.106.26.042	42
11.106.26.044	44

PRODUCT CODE	LENGTH
11.106.26.046	46
11.106.26.048	48
11.106.26.050	50
11.106.26.052	52
11.106.26.054	54
11.106.26.056	56
11.106.26.058	58
11.106.26.060	60
11.106.26.062	62
11.106.26.064	64
11.106.26.066	66

PRODUCT CODE	LENGTH
11.106.26.068	68
11.106.26.070	70
11.106.26.072	72
11.106.26.074	74
11.106.26.076	76
11.106.26.078	78
11.106.26.080	80
11.106.26.082	82
11.106.26.084	84
11.106.26.086	86
11.106.26.088	88

Tibia Cancellous Locking Screw

- Full Thread
- Self Tapping
-  SW4.8



Ti:

PRODUCT CODE	LENGTH
11.106.27.025	25
11.106.27.030	30
11.106.27.035	35
11.106.27.040	40
11.106.27.045	45

PRODUCT CODE	LENGTH
11.106.27.050	50
11.106.27.055	55
11.106.27.060	60
11.106.27.065	65
11.106.27.070	70

PRODUCT CODE	LENGTH
11.106.27.075	75
11.106.27.080	80
11.106.27.085	85
11.106.27.090	90

End Cap

-  SW3.5

TI:

PRODUCT CODE	TYPE	LENGTH
11.106.29.000	STANDARD-LONG	18
11.106.29.008	STANDARD	12
11.106.28.005	EXTRA-LONG	18
11.106.28.010	EXTRA-LONG	22
11.106.28.015	EXTRA-LONG	28

